

# NCLwipes

## Disinfectant Wipes

Lemon Fresh

Wipes



NCL Product Technical Data

### DESCRIPTION

NCLwipes™ Disinfectant Wipes clean, disinfect and deodorize in one easy step. NCLwipes™ Disinfectant Wipes save time and labor and are easy to use. No mixing, No measuring, No diluting. Simply pull a wiper from the canister, wipe the surface and discard.

NCLwipes™ Disinfectant Wipes is a broad spectrum disinfectant that kills HIV-1 (AIDS virus) in 1 minute, Norwalk Virus in 10 minutes, and most other viruses, bacteria in 4 minutes.

### PRODUCT FEATURES / BENEFITS

- Kills 99.9% of Bacteria in just 15 seconds.
- Kills Cold and Flu Virus\*<sup>1</sup>.
- Kills Pseudomonas, Staph, Salmonella, MRSA, CA-MRSA, E.coli, VRE, Norovirus, HBV, HCV in 4 minutes!
- Cleans, Disinfects and Deodorizes.
- Fast 4 minute disinfectant contact time
- Broad spectrum kill claims of bacteria and viruses
- "One-Step" disinfecting on hard, non-food, non-porous surfaces.
- Controls the growth of mold, mildew and odors.
- Solvent free, bleach free, alcohol free.
- Quat-Based, Alternative to alcohol or bleach-based wipes.
- Non-flammable
- Lemon Fresh Fragrance
- Strong, Absorbent Wiper does an excellent job of cleaning
- Cloth like texture for scrubbing power
- Sturdy, Attractive, Functional Point-of-use Canister
- Factory Sealed, Tamper proof
- Pull, wipe, discard in trash.
- Recyclable packaging
- Streak free.
- Reduces cross contamination on hard, non-porous treated surfaces
- Designed for the toughest tasks.
- Ready to use, does not require mixing or diluting
- No hunting for spray bottles or towels.

### APPLICATION AREAS

Schools, Colleges, Hospitals, Medical Offices, Nursing homes, Day care centers, Nurseries, Institutional facilities, USDA inspected food processing facilities, Office buildings, Restaurants, Grocery stores, Convenience stores, Hotels, Motels, Prisons, Salons

### APPLICATION SURFACES

Chairs, Counters, Desks, Doorknobs, Faucets, Glass, Light switches, Showers, Sinks, Tables, Telephones, Toilets, Diaper changing tables, and other hard nonporous surfaces made of acrylic, metal, plastic, stainless steel, upholstery- vinyl and plastic.

### DIRECTIONS FOR USE

**TO CLEAN AND DEODORIZE:** Wipe surface with wipe. Let air dry. Discard used wipe in trash. Do not flush in toilet.

**TO CLEAN AND DISINFECT:** Wipe surface with wipe until surface is visibly wet. Use enough wipes to thoroughly wet surfaces. Use enough wipes to keep surfaces visibly wet for 4 minutes. Allow surface to remain wet for 4 minutes. For *Candida albicans*, treated surfaces must remain wet for 8 minutes. For *Norovirus* and *Rotavirus*, treated surfaces must remain wet for 10 minutes. If surfaces are extremely dirty, clean first with another wipe before disinfecting. Discard used wipe in trash. Do not flush in toilet.

Rinse food contact surfaces with clean, potable water after disinfecting. Do not use to disinfect dishes, glassware or utensils.

#### FOR MOLD AND MILDEW

NCLwipes™ Disinfectant Wipes will effectively inhibit the growth of mold and mildew and the odors caused by them when applied to hard, nonporous surfaces. Follow disinfection instructions. Repeat treatment every seven days, or more often if new growth appears.

#### TO SANITIZE

To sanitize hard, non-porous, non-food contact surfaces; wipe surface, using enough wipes to allow treated surface to remain wet for 15 seconds for *Staphylococcus aureus*, *Klebsiella pneumoniae*, and *Campylobacter jejuni*. Let surface dry. For highly soiled surfaces clean excess dirt first.

**\*KILLS HIV, HBV AND HCV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), or Hepatitis B Virus (HBV) or Hepatitis C Virus (HCV).

#### SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 or HBV or HCV ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS.

**Personal Protection:** Cleanup must always be done wearing protective gloves, gowns, masks and eye protection.

**Cleaning procedure:** Blood and other body fluids containing HIV, HBV or HCV must be thoroughly cleaned from surfaces and objects before application of NCL® Wipes Disinfectant Wipes. NCLwipes™ Disinfectant Wipes can be used for this purpose.

**Contact time:** Leave surface wet for 1 minute for HIV-1. Use appropriate contact time as listed on label for all other viruses, bacteria and fungi claimed.

**Disposal of Infectious Materials:** Blood and other body fluids, cleaning materials and clothing must be disposed of according to local regulations for infectious waste disposal.

This product is for use on hard nonporous surfaces not conducive to treatment by immersion or excess liquid against Avian Influenza A.

\*<sup>1</sup> On hard, non-porous, non-food surfaces.

continued on back side



World Class Cleaning & Hygiene Solutions™

TDS-4542-0120



# NCL *wipes*

## Disinfectant Wipes

*Lemon Fresh*

Wipes

### EFFICACY

Bactericidal Test Results	Contact Time (in min.)
Acinetobacter baumannii	4
Brevibacterium ammoniagenes	4
Campylobacter jejuni	4
Enterobacter aerogenes	4
Enterococcus faecalis	4
Enterococcus faecalis – Vancomycin resistant	4
Escherichia coli	4
Escherichia coli O157:H7	4
Escherichia coli ESBL (antibiotic resistant)	4
Klebsiella pneumoniae	4
Legionella pneumophila	4
Listeria monocytogenes	4
Pseudomonas aeruginosa	4
Pseudomonas cepacia – Strain 25416	4
Salmonella enterica	4
Salmonella schottmuelleri	4
Salmonella typhi	4
Serratia marcescens	4
Shigella dysenteriae	4
Staphylococcus aureus	4
Staphylococcus aureus – Multi-drug resistant	4
Staphylococcus aureus – MRSA	4
Staphylococcus aureus – CA-MRSA (USA 300/400)	4
Staphylococcus aureus – VISA	4
Streptococcus pyogenes	4
Vibrio cholerae	4
Virucidal Test Results	Contact Time (in min.)
*Hepatitis B Virus (HBV)	4
*Hepatitis C Virus (HCV)	4
*Herpes Simplex Virus Type 1	4
*Herpes Simplex Virus Type 2	4
*HIV-1 (AIDS Virus)	1
*Human Coronavirus	4
*Influenza Type A Virus / Brazil (the Flu Virus)	4
*Influenza A Virus H1N1	4
*2013 Influenza A virus (H7N9)	4
*Norwalk Virus (Norovirus)	10
*Respiratory Syncytial Virus	4
*Rotavirus	10
*SARS associated Coronavirus	4
*Vaccinia Virus	4
Virucidal Test Results (Animal)	Contact Time (in min.)
Avian Influenza Virus (H3N2)	4
Avian Influenza Virus (H5N1)	4
Canine Distemper Virus	4
Feline Calicivirus	10
Newcastle's Disease Virus	4
Pseudorabies Virus	4
Fungicidal Test Results	Contact Time (in min.)
Candida albicans	8

Respiratory illnesses attributable to Pandemic 2009 H1N1 are caused by influenza A virus. This product is a broad-spectrum hard surface disinfectant that has been shown to be effective against influenza A virus tested and listed on the label and is expected to inactivate all influenza A viruses including Pandemic 2009 H1N1 (formerly called swine flu).

### SPECIFICATIONS

Wipe Size .....	7" x 8"
Wipes per Canister .....	125
Canisters per Case .....	6
Net Wt / Canister .....	2.1 lbs
Color .....	White Cloth
Fragrance .....	Lemon Fresh
pH .....	7.0 ± 0.5
Flammability .....	Not Flammable
Shelf Life .....	1 year minimum in original unopened container

### SAFETY INFORMATION

#### PRECAUTIONARY STATEMENTS

#### HAZARD TO HUMANS AND DOMESTIC ANIMALS

**Caution:** Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet.

#### FIRST AID:

If in eyes: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after first 5 minutes, then continue rinsing eye.

Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling poison control center or doctor or going for treatment.

\*1 On hard, non-porous, non-food surfaces.



World Class Cleaning & Hygiene Solutions™

401 N. 10th Street • Philadelphia, PA 19123 USA

Toll Free: (800) 628-2436

Phone: (215) 922-1200 • Fax: 215-922-5517

Website: [www.nclonline.com](http://www.nclonline.com)

Contact E-mail: [info@nclonline.com](mailto:info@nclonline.com)