

NCLwipes

Disinfectant Wipes

Lemon Fresh

Wipes



NCL Product Technical Data

DESCRIPTION

NCLwipes™ Disinfectant Wipes clean, disinfect and deodorize in one easy step. **NCLwipes™ Disinfectant Wipes** save time and labor and are easy to use. No mixing, No measuring, No diluting. Simply pull a wipe from the canister, wipe the surface and discard. **NCLwipes™ Disinfectant Wipes** is a broad spectrum disinfectant that is effective against SARS-Related Coronavirus 2¹ (SARS-CoV-2) in 15 seconds, Norwalk Virus in 10 minutes, and most other viruses, bacteria in 4 minutes.

PRODUCT FEATURES / BENEFITS

- Effective against SARS-Related Coronavirus 2¹ (SARS-CoV-2) in 15 seconds.
- Kills Cold and Flu Virus.^{1,1}
- Controls the growth of odor-causing bacteria for 72 hours.
- 3 steps in 1- wipes, scrubs stubborn spots and disinfects² so you finish faster.
- Everyday surface¹ disinfectant.
- Kills 99.9% of Bacteria in just 15 seconds.³
- Effective in the presence of 5% serum contamination.
- Eliminates 99.9% of Pseudomonas aeruginosa, Staphylococcus aureus - Methicillin-Resistant, Escherichia coli, Enterococcus faecalis- Vancomycin Resistant, Salmonella enterica in just 4 minutes!
- Leaves surfaces¹ clean and disinfected.
- An effective way to clean and deodorize against odor-causing bacteria.
- Faster acting wipe formula.
- A one-step disinfectant² that is effective against a broad spectrum of bacteria when used as directed.
- Phosphate free.
- This product will clean and disinfect hard, non-porous, non-food surfaces.
- Lemon Fresh Fragrance.
- Strong wipe for tougher jobs.
- Textured for better scrubbing power.
- Just pull wipe and toss.
- Recyclable packaging.
- Streak free.
- Reduce risk of cross contamination between treated surfaces.¹
- Ready to Use.
- Handles daily clean-ups with no need for additional sprays, sponges or buckets.

APPLICATION AREAS

Schools, Colleges, Hospitals, Medical Offices, Nursing homes, Day care centers, Nurseries, Institutional facilities, USDA inspected food processing facilities, Office buildings, Restaurants, Grocery stores, Convenience stores, Hotels, Motels, Prisons, Salons

APPLICATION SURFACES

Chairs, Counters, Desks, Doorknobs, Faucets, Glass, Light switches, Showers, Sinks, Tables, Telephones, Toilets, Diaper changing tables, and other hard nonporous surfaces made of Acrylic, Metal, Plastic, Stainless steel, Upholstery-vinyl and Plastic.

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

Do not use as a diaper wipe or for personal cleansing. This is not a baby wipe! Do not flush in toilet.

TO CLEAN AND DEODORIZE: Wipe hard, nonporous surface with wipe. Let air dry. Discard used wipe in trash. Do not flush in toilet.

TO CLEAN AND DISINFECT: Wipe hard, nonporous surface with wipe until surface is visibly wet. Use enough wipes to visibly wet surfaces. Use enough wipes to keep surfaces visibly wet for 4 minutes. Allow surface to remain visibly wet for 4 minutes. **For Norwalk Virus, Feline Calicivirus, Norovirus and Rotavirus**, treated surfaces must remain visibly wet for 10 minutes. If surfaces are visibly dirty, clean first with another wipe before disinfecting. Discard used wipe in trash. Do not flush in toilet.

FOR MOLD AND MILDEW STAIN CONTROL

NCLwipes™ Disinfectant Wipes will effectively inhibit the growth of mold and mildew and the odors caused by them when applied to hard, nonporous surfaces. Wipe surface and let it remain visibly wet for 4 minutes. Repeat treatment every seven days, or more often if new growth appears.

SANITIZING

To sanitize hard, non-porous, non-food contact surfaces; wipe surface, using enough wipes to allow treated surface to remain visibly wet for 15 seconds for Staphylococcus aureus, Klebsiella pneumoniae, and Campylobacter jejuni. Let surface dry. For visibly soiled surfaces clean visible soil first.

***KILLS HIV, AND HBV AND HCV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), or Hepatitis B Virus (HBV) or Hepatitis C Virus (HCV).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 or HBV or HCV ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS.

Personal Protection: Clean-up must always be done wearing protective gloves, gowns, masks and eye protection.

Cleaning procedure: Blood and other body fluids containing HIV or HBV or HCV must be visibly cleaned from surfaces and objects before application of **NCLwipes™ Disinfectant Wipes**. **NCLwipes™ Disinfectant Wipes** can be used for this purpose.

Contact time: Leave surface visibly wet for 1 minute for HIV-1. Use appropriate contact time as listed on label for all other viruses and bacteria claimed.

Disposal of Infectious Material: Blood, body fluids, cleaning materials and clothing must be disposed of according to local regulations for infectious waste disposal.

This product is for use on hard nonporous surfaces not conducive to treatment by immersion or excess liquid against Avian Influenza A.

This product is not for use on medical device surfaces.

*1 On hard, non-porous, non-food surfaces. ² When used according to the disinfectant directions ³ When used according to the sanitization directions

continued on back side



World Class Cleaning & Hygiene Solutions™

TDS-4542-0124

